

Christina Do, DDS
Costa Mesa, CA
glidewell.io user since 2018



Michael Gerlach, DDS
Conroe, TX
glidewell.io user since 2019



Steven Bui, DDS
Westminster, CA
glidewell.io user since 2020



Damanjeet Singh, DDS
West Palm Beach, FL
glidewell.io user since 2019



SAME-VISIT RESTORATIONS FOR EVERY DENTIST

FEATURING BRUXZIR® NOW ZIRCONIA



glidewell.io™
IN-OFFICE SOLUTION

www.glidewell.io ■ 844-949-7184

SAME-VISIT BRUXZIR® CROWNS ARE NOW A REALITY



With [glidewell.io™](#), you can quickly design and mill crowns from the most clinically proven zirconia in dentistry.

The [glidewell.io™](#) In-Office Solution is the only way to mill BruxZir® crowns right in your practice — in fewer steps and less time than other systems. Two out of three patients would change dentists for same-visit crowns.¹ Set your practice apart with single-appointment, lab-quality restorations, with manufacturer-direct support from the largest dental laboratory in the world.

[glidewell.io](#) is available at half the price of competitor systems — and offers twice the output — giving your practice an edge in same-visit dentistry.



BruxZir® NOW

The zirconia you trust, ready to deliver straight from the mill — just cement and move on.

Mills in under 40 minutes. \$35 per crown.*

[glidewell.io](#) is the fastest-growing chairside CAD/CAM system. Hear your peers share why.



BruxZir® NOW Zirconia Crowns

Boost your productivity with the only in-office zirconia that lets you skip the oven wait.

CrownAI™ One-Click Design

Let the software do the work — AI-powered designs mean you don't have to be a tech.

1. Exevia GmbH – German dental patient survey, November 2015.

* Milling time based on customer usage data. Gold subscription service pricing.

ONE APPOINTMENT. FOUR SIMPLE STEPS.

Your Laboratory Partner Is Now Chairside

The glidewell.io™ In-Office Solution encompasses the flexibility, simplicity and efficiency clinicians look for when considering a CAD/CAM system. The No Oven Wait (NOW) materials exclusive to glidewell.io are all ready to place straight out of the mill, skipping the oven time required for other in-office solutions.

1. SCAN



The glidewell.io process begins with a digital impression taken with the intraoral scanner of your choice.

2. DESIGN



Next, you simply confirm the auto-generated crown design created by proprietary AI technology.

3. MILL



Once you've approved your design and loaded the block, milling begins.

4. DELIVER



You're set for delivery, all in under one hour!

Real-Time, Direct Support from Glidewell

With glidewell.io, you get one-call resolution from an expert technical support team that can help you with your restoration design and remote into your system to make adjustments. Our technicians are available via text, chat, email or phone.

FAST, ACCURATE RESULTS POWERED BY AI



fastdesign.io™ at a Glance

The fastdesign.io™ Software and Design Station provides you with ready-to-use designs auto-generated by Glidewell's CrownAI™ technology. All you have to do is complete your margins and approve the design.

- Computer-assisted margin marking based on intraoral scan
- Highly accurate, individualized designs
- Ready to mill in less than 5 minutes



Watch Demo Video



The designs that come back from the software are very precise and don't need much, if any, adjustment. It's all very straightforward and fast.

- Michael Buck, DDS | Creston, IA



The Superiority of CrownAI™


Glidewell's CrownAI™ uses exclusive artificial intelligence (AI) and machine-learning algorithms to ensure the highest accuracy. Pulling data from Glidewell's extensive database of past cases, the software automatically combines that information with the case at hand and provides you with a 3D design ready for your approval. No more designing from scratch; no more inaccurate library teeth.




margin marked



AI design



milling preview



AI does the work — so you can spend your time with your patients.

SMALL FOOTPRINT, BIG VALUE AND BRUXZIR® CROWNS



fastmill.io™ at a Glance

The fastmill.io™ In-Office Mill is smaller than most microwave ovens, but it can deliver big-time results for your practice. It's the only chairside milling unit capable of producing BruxZir® NOW crowns.

- Mills fully sintered crowns; no oven required
- Plug-and-play operability
- Manufactured in the USA by Glidewell



Watch Demo Video



glidewell.io is the simplest, most well-thought-out CAD/CAM system yet – great for practices large and small.

- Brian Wilkes, DDS | Chalfont, PA



Fast and Efficient Dentistry

The power behind the fastmill.io eliminates the need for an oven wait, reducing appointment times and increasing your productivity. While the crown is milling for one patient, you have the flexibility to move on to your next patient, allowing for a more efficient workflow. Just leave the mill running and go about the rest of your routine.

In under 40 minutes, your BruxZir crown will be milled and ready for delivery when you return to the other patient.

Scheduling Example



CHAIR 1

- 1. Prep: Doctor
- 2. Scan: Assistant
- 3. Design: Assistant
- 4. Mill: 40 Minutes
- 5. Polish: Assistant
- 6. Delivery: Doctor


BruxZir NOW
crowns ready to deliver in under 45 minutes

CHAIR 2


- 1. Prep: Doctor
- 2. Scan: Assistant
- 3. Design: Assistant
- 4. Mill: 40 Minutes
- 5. Polish: Assistant
- 6. Delivery: Doctor

3

Things That Are Bigger Than the fastmill.io



microwave



laser printer



every other mill on the market

THE ZIRCONIA YOU TRUST. READY TO DELIVER.



BruxZir® NOW

**BruxZir®
esthetic NOW**

BruxZir® NOW at a Glance

Exclusive to Glidewell milling systems, BruxZir® NOW offers unmatched strength and efficiency. By eliminating the need for oven time and post-processing, BruxZir NOW crowns are ready to go right from the mill.

- Saves you time and money by eliminating the oven time required for other in-office zirconias
- Also available in BruxZir Esthetic formulation for anterior cases
- Ideal for posterior crowns demanding maximum strength



Watch BruxZir NOW
Case Video



A lot of dentists say, 'I'm not going to bother with chairside milling if I can't do BruxZir restorations.' Well, with glidewell.io, you can have it both ways.

- Scott Barton, DDS | San Antonio, TX



Save Time and Increase Productivity with BruxZir NOW

BruxZir NOW Zirconia gives you unrivaled strength while reducing the steps required to mill same-visit crowns. Unlike other in-office zirconias, BruxZir NOW can be delivered straight from the mill, saving you valuable time.

BRUXZIR® NOW ZIRCONIA CHAIRSIDE WORKFLOW



Comparable in-office materials are not fully sintered, which means they must be baked after milling, requiring the doctor to interrupt their workflow multiple times for post-processing steps. Instead, with the NOW materials exclusive to glidewell.io, you can just let the mill do its job and simply come back when the crown is ready to deliver — no oven wait required.

CHAIRSIDE WORKFLOW FOR OTHER IN-OFFICE ZIRCONIAS



BruxZir® NOW



BEFORE



AFTER

Failed restorations due to recurrent caries were replaced with high-strength BruxZir® NOW crowns that will last long into the future.

BruxZir® esthetic NOW



BEFORE



AFTER

New BruxZir Esthetic NOW restorations replaced these failed PFM crowns and brought lifelike esthetics back to the patient's smile.


IDEAL RESTORATIVE MATERIALS FOR EVERY INDICATION




BruxZir® NOW



BruxZir®
esthetic NOW



CAMouflage® NOW



BioTemps® NOW



Obsidian®

Material Comparison Chart

		BruxZir® NOW	BruxZir® Esthetic NOW	CAMouflage® NOW	BioTemps® NOW	Obsidian® Milling Blocks
Material Composition		Solid zirconia	Solid zirconia	Ceramic-filled nanohybrid resin composite	PMMA — Poly(methyl methacrylate)	Monolithic lithium silicate ceramic
Milling Time		< 40 minutes	< 45 minutes	< 10 minutes	< 10 minutes	< 20 minutes
Oven Needed						✓
Ideal Uses	Anterior		✓	✓	✓	✓
	Posterior	✓		✓	✓	✓
	Inlay/Onlay			✓		✓
	Veneers			✓		✓
	Bridges	Coming soon!			✓	
	Temporary Solution				✓	



With glidewell.io™, we're providing patients with the best materials — all with the convenience of single-visit dentistry.

- Eric Seidel, DMD | Gettysburg, PA



BioTemps® NOW



BEFORE



AFTER

BioTemps® NOW provisional restorations served as well-fitting temporaries for #12-15 while the patient awaited final treatment.

BioTemps® NOW BRIDGE BLOCK



BEFORE



AFTER

This BioTemps NOW bridge was milled in-office and functioned well as a temporary solution while protecting the patient's teeth during healing.

CAMouflage® NOW



BEFORE



AFTER

CAMouflage NOW overlays on teeth #30-31 remedied the failing composite restorations on those teeth.

Obsidian® MILLING BLOCKS




BEFORE



AFTER

Single-visit Obsidian crowns on teeth #7-8 replaced failing PFMs, and single-visit Obsidian veneers improved the symmetry on #9-10.

OPTIMIZE YOUR SCANS. STREAMLINE YOUR WORKFLOW.



iTero Element® 5D Plus at a Glance

The iTero Element® 5D Plus optimizes your intraoral scans while opening the door to the most efficient digital workflow Glidewell offers.

- Experience 3D imaging, intraoral color and near-infrared (NIRI) technology
- Revolutionize long-term treatment with TimeLapse technology
- Encourage patient engagement with clear, easy-to-understand visuals

AN EASY ENTRY POINT — WHATEVER YOUR SCANNER



*glidewell.io is a perfect way to begin milling in-office restorations.
Plus, it integrates easily with our intraoral scanner.*

- Susan McMahon, DMD | Pittsburgh, PA



CEREC®
Primescan



3shape
TRIOS



Carestream
DENTAL



MEDIT i700



3M True
Definition

Open System at a Glance

With the glidewell.io™ In-Office Solution, you can practice same-visit dentistry your way. While glidewell.io integrates seamlessly with the iTero Element® line, it also works with every popular scanner on the market today. So if you already own an intraoral scanner, or prefer a different brand, we can help you link it with your glidewell.io system and achieve an efficient workflow.

Using a CEREC mill? 1 out of 5 glidewell.io customers are CEREC trade-ins.

Enjoy effortless connectivity with:

- CEREC® Primescan
- 3Shape TRIOS®
- Carestream™ CS 3600
- Medit i700
- 3M™ True Definition
- and many more!

QUESTIONS? HEAR WHAT USERS HAVE TO SAY



What if I need help with one of my glidewell.io™ cases?

“We can always reach out to Glidewell support for help if needed. The support team is easy to access, too. A member of our staff just makes a quick phone call and the Glidewell tech can jump right into our system and adjust a couple of things in our design. Then we can load up a block and start milling just a few minutes later.”

– Thomas Rolfes, DDS
Costa Mesa, CA

How can my team and I learn to master our glidewell.io system?

“Glidewell never stops. When I first purchased glidewell.io, they gave us onsite training and resources for how to get started. Since then, they’ve never stopped giving me things that make my dentistry better and more exciting.”

– Thomas Fuchs, DMD | Louisville, KY



Would investing in glidewell.io be a wise financial move for my practice?

“glidewell.io has boosted production at our office so much that we’ve purchased two additional mills to keep up with the demand for same-visit crowns. Each mill has added approximately \$1,000 in revenue to our practice per day. Investing in glidewell.io is the best financial decision we’ve ever made.”

– Joshua Prentice, DDS | Bend, OR



Calculate Your ROI

After just four months of ownership, offices using glidewell.io report a 65% or more decrease in monthly lab-related expenses.* Scan the code to calculate your ROI or visit www.glidewell.io/roi-calculator.

*Based on Glidewell internal data.

WORLD-CLASS SUPPORT FROM RESTORATIVE EXPERTS

Direct Support from Glidewell

Unlike systems sold by a distributor, glidewell.io puts your dental lab chairside, so you get maximum value from your investment.

- Hands-on installation and one-on-one training in your practice on Day One
- Immediate access to glidewell.io support through text, chat, email or phone
- Real-time remote design assistance from skilled support technicians



Continuing Education

If you're ready for a deeper understanding of restoration design and more complex cases, glidewell.io provides learning opportunities to guide the way.

- Live, hands-on advanced training courses
- Clinician-led monthly webinars
- FREE on-demand educational articles and videos

DON'T GET LEFT BEHIND

A BruxZir®
NOW crown
is milled every

3.75

minutes by a
dentist somewhere
near you

*Percy Torkornoo, DDS | Fairborn, OH
glidewell.io user since 2019*



*Lovedeep Randhawa, DMD | Surrey, British Columbia
glidewell.io user since 2019*



*Davis Gardner, DDS | Chesterfield, VA
glidewell.io user since 2018*



Ready to learn more about what glidewell.io
can add to your practice?

Let's talk! Scan the code and provide your information and we'll be in touch soon.

glidewell.io™
IN-OFFICE SOLUTION

www.glidewell.io | 844-949-7184